

## **Rigidscan LR 7/LR 7s** **Rigid Endoscope Tester**

**An easy to use  
device for visual  
inspection of rigid  
endoscopes**



**The LR 7 is the optimal instrument  
for quality checking in the processing  
of instruments as well as for evaluation purposes**

### **NEW!**

**There is now a special version available for checking short  
endoscopes - the LR7s!**

**The following optical damage may be assessed:**

- broken lenses
- dirt on lenses
- bubbles and discolouration in cement

### **Range of application**

#### **processing / sterilization**

- check for damage
- confirm image quality

#### **medical technician, checking for**

- quality on delivery
- damage assessment
- defining origin of malfunction

## Specifications\*

diameter of endoscopes:	suitable for all sizes (different adapters available)
eyepieces:	eyepieces according to DIN 58105
length of endoscopes:	Universal
endoscope types:	all well-established types
calibration:	not necessary

### **optical properties**

form of lens:	biconvex
magnification:	25 dpt
Lens diamètre:	19 mm

**\*Specifications are subject to change without notice**